



Sodium Periodate Powder Pillows, 10 mL, pk/100

Ürün Numarası:

2107769

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



Powder reagent in sealed packages for the determination of Manganese, cold periodate method. For 10 mL sample.

For determination of manganese by the periodate oxidation method. EPA approved for reporting wastewater analyses (digestion required). Hach method 8034. Sample size: 10 mL Pack of 100 Powder Pillows.

Teknik özellikler

Hacim / Paket boyutu:	100 /pk
Metot:	Periodate Oxidation
Platform:	Powder Pillows
Test sayısı:	100