



Be Right™



Sodium periodate powder pillows for manganese, 25 mL, pk/100

Ürün Numarası:

98499

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



Powder reagent in sealed package for the determination of Manganese, cold periodate method, and Iron titrations. For 25 ml sample.

For determination of manganese by the periodate oxidation method. EPA approved for reporting wastewater analyses (digestion required). Hach method 8034. Sample size: 25 mL Pack of 100 Powder Pillows.

Teknik özellikler

Hacim / Paket boyutu: 100 /pk

Metot: Periodate Oxidation

Platform: Powder Pillows

Test sayısı: 100