



Be Right™



Digital titrator cartridge, hydrochloric acid , 8.00 N

Ürün Numarası:

1439001

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



For base concentration determination by Digital Titrator titration. Digits x 0.01 = meq

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Teknik özellikler

Metot:	Phenolphthalein End Point
Metot numarası:	8233 (Base)
Ölçüm aralığı:	10 - 4000 meq/L Base
Parametre:	Acid-Base
Platform:	Dijital Titratör
Test sayısı:	40 - 100
Titrant:	8.00 N HCl

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)