




## CDTA Digital titrator cartridge, 0.0800 M

|                |             |
|----------------|-------------|
| Ürün Numarası: | 1440201     |
| TRY Fiyat:     | Bize Ulaşın |
| Mevcut         |             |

Cyclohexanediameine tetraacetic acid standard solution. For Calcium and Total Hardness determination by Digital Titrator titration. 0.0800 M. Digits x 0.1 = mg/L CaCO<sub>3</sub> in 100 mL. Digits x 0.0001 = mmole.

### Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

### Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

### Versatile

More than 50 parameter/range combinations are available

### Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

---

## Teknik özellikler

|                |                                 |
|----------------|---------------------------------|
| Metot:         | CDTA Titration                  |
| Ölçüm aralığı: | 10 - 160 mg/L CaCO <sub>3</sub> |
| Platform:      | Dijital Titrator                |