



## Reagent, hydrazine, AccuVac ampules, pk/25

Ürün Numarası:

2524025

Tehlikeli

**ÜRETİMDEN KALDIRILAN ÜRÜN**  
Bu ürün artık mevcut değil.



Önerilen yedek parçalar:

- 179032

For determination of Hydrazine by the p-Dimethylaminobenzaldehyde method, using AccuVac Ampules. Hach Method 8141. Range: 4 - 600 µg/L. Pack of 25 ampules.

### Fast

AccuVac Ampules provide quick results. Just snap, fill and read - saving you valuable time.

### Easy

No need to measure sample or reagents. The evacuated ampules draw in the correct sample volume, so anyone can obtain reliable results with minimal effort.

### Accurate

Optical quality evacuated vials assure precise measurements and reagent stability.

### Use them anywhere you need results.

Ampules are both reagent containers and photometric measurement cells. No need to transport reagents or samples.

## Teknik özellikler

Cihaz:	All except DR/850, DR/820, PC II
Metot:	p-Dimethylaminobenzaldehyde
Metot numarası:	8141
Ölçüm aralığı:	4 - 600 µg/L N <sub>2</sub> H <sub>4</sub>
Parametre:	Hydrazine
Platform:	Accuvac
Test sayısı:	25

## Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)