



Acid reagent powder pillows, 25 mL, pk/1000

Ürün Numarası:
TRY Fiyat:
Sevk tarihi gösterilmedi

104228
Bize Ulaşın

Tehlikeli



For determination of high range silica by the silicomolybdate method (Hach method 8185) Range: 1-100.0 mg/L SiO₂. Also used in the iodometric method (Hach method 10100) for the determination of high range hypchlorite (50-150 g/L as Cl₂). Pack of 1000 powder pillows.

Teknik özellikler

Metot: Silicomolybdate
Platform: Powder Pillows

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)