



**Be Right™**



## Test strips free chlorine, 0-600 mg/L, 100 tests

Ürün Numarası:

2890200

TRY Fiyat:

Bize Ulaşın

Mevcut

Hach test strips are manufactured with a unique reagent pad at the tip, using proven chemistries based on standard reference methods. Just dip the strip, wait until color development occurs, and compare the color of the reacted pad to the chart printed on the bottle.

**Ideal for anyone performing semi-quantitative pre-screenings, spot-checks, or field tests**

**Simplify testing and get results in less than 1 minute, in the field or lab**

**Eliminate chemical handling and clean-up**

**Reduce guesswork by screening samples**

**\*Not for use in sea or brackish water**

### Teknik özellikler

Birim:	100 tests
Dipnot:	*mg/L unless otherwise noted; ppb = µg/L; ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
En küçük artış adımları:	Steps: 0, 25, 50, 100, 200, 400, 600
Ölçüm aralığı:	0 - 600 mg/L Cl <sub>2</sub>
Parametre:	Free Chlorine, High Range
Platform:	Test Strip