



Filter, glass fiber, diameter 55 mm, preweighed standard methods TSS procedure, 100/pk

Ürün Numarası:

2546120

TRY Fiyat:
Mevcut

Bize Ulaşın

Individually preweighed filter papers, preweighed in compliance with Standard Methods Total Suspended Solids procedure. Each in a planchet with tare weight listed. Each tare weight meets Standard Methods specifications.

Teknik özellikler

Birim:	100/pk
Çap:	55 mm
Gözenek boyutu:	1.5 µm
Hacim / Paket boyutu:	100 /pk
Malzeme:	Glass Microfiber