



MEL/P-A safe drinking water lab

Ürün Numarası:

2569600

TRY Fiyat:

Bize Ulaşın

Teslimat tarihi için arayın

Tehlikeli



MEL/P-A Safe Drinking Water Lab

Economical portable lab monitors key water quality parameters plus Total Coliforms and E. Coli (Presence-Absence Method). MEL P/A Safe Drinking Water Lab includes portable incubator, reagent sets, portable long-wave UV lamp, manual and apparatus.

Consumables only

Choose media and order seperately

Nitrate and Chlorine determined using visual color discs

Portable incubator used to determine total and fecal coliforms directly in the field

Teknik özellikler

Hassasiyet:

CFU/mL