



pH test kit, 5.5 - 8.5 pH, model 17F

Ürün Numarası:

147006

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



Simple, low-cost color disc method is the most accurate visual test kit method

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Teknik özellikler

Metot:	Color disc/ Bromthymol blue
Model:	17F
Muhafaza Türü:	B
Ölçüm aralığı:	5,6 - 8,4 pH
Parametre:	pH, mid range
Platform:	Renk Diski
Smallest increment:	0.1 pH
Test kiti modeli:	17F
Test sayısı:	200