



Be Right™



Test kit, phosphorus, orthophosphate (reactive) , model PO-19A

Ürün Numarası:

224801

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



Simple, low-cost color disc method is the most accurate visual test kit method

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Teknik özellikler

Metot:	Color disc/Ascorbic acid
Model:	PO-19A
Muhafaza Türü:	F
Ölçüm aralığı:	Triple Range
Ölçüm aralığı (2):	0 - 4 mg/L PO ₄
Ölçüm aralığı (3):	0 - 40 mg/L PO ₄
Parameter:	Phosphorus
Parametre:	Phosphorus, Orthophosphate (Reactive) - as PO ₄ ³⁻
Platform:	Renk Diski
Smallest increment:	0.02
Standart aralık:	0 - 0.8 mg/L PO ₄ ; 0 - 4 mg/L PO ₄ ; 0 - 40 mg/L PO ₄
Test kiti modeli:	PO-19A
Test sayısı:	(Includes filtration for turbid samples.)
100	