



Chloride Low Range Test Kit, Model 8-P

Ürün Numarası:

144001

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



Chloride, single parameter test kit, drop count titration (range: 5 - 100

mg/L) / silver nitrate (range: 20 - 400 mg/L), blue polypropylene case

Using silver nitrate titrant, this kit provides fast, accurate detection of chloride.

Low cost and easy to use drop count titration method using a silver nitrate titrant

Ideal for quickly determining quantitative results

Teknik özellikler

Metot:	Drop count titration/
	Silver nitrate
Muhafaza Türü:	B
Ölçüm aralığı:	5 - 100 mg/L & 20 - 400 mg/L Cl ⁻
Parametre:	Chloride - As Cl ⁻
Platform:	Drop Count Titration
Smallest increment:	5
Test kiti modeli:	8-P
Test sayısı:	100