



## Hardness (Total & Calcium) Test Kit, Model HA-4P MG-L

Ürün Numarası:

145701

TRY Fiyat:  
Mevcut

Bize Ulaşın

Tehlikeli



Hardness, Total & Calcium (Magnesium by difference) drop count titration/EDTA test kit, range 20 - 400 mg/L CaCO<sub>3</sub>, 100 tests.

The Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field.

**Utilises drop count titration - a simple, economical titration method**

### Teknik özellikler

Kutudakiler:	Includes five reagents, bottle, measuring tube, instruction sheet, and carrying case.
Metot:	Drop count titration/EDTA
Muhafaza Türü:	D
Ölçüm aralığı:	20 - 400 mg/L CaCO <sub>3</sub>
Parametre:	Hardness, total & Calcium - (Magnesium by difference)
Platform:	Drop Count Titration
Smallest increment:	20
Test kiti modeli:	HA-4P
	MG-L
Test sayısı:	100

### Kutudakiler

Includes five reagents, bottle, measuring tube, instruction sheet, and carrying case.