



Chlorine (total) test kit, model CN-21P, 10 - 200 mg/L, 100 tests

Ürün Numarası:

2444400

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



Method / Chemistry: Drop count titration / Thiosulfate

Chlorine Test Kits are manufactured by Hach Company to ensure compliance with stringent quality control standards. Using the Drop Count Titration method, the CN-21P Chlorine Test Kit is used to determine high concentrations of chlorine in sanitizing applications. 10 - 200 mg/L; contains 100 tests.

Simple, economical titration method

Kit contains all required reagents, solutions, and glassware

Teknik özellikler

Dipnot:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Metot:	Drop count titration/ Thiosulfate
Muhafaza Türü:	B
Ölçüm aralığı:	10 - 200 mg/L
Parametre:	Chlorine, mid and high range - As total Cl ₂
Platform:	Drop Count Titration
Smallest increment:	10
Test kiti modeli:	CN-21P
Test sayısı:	100