



Test kit, total chlorine, model CN-DT, 20-2000 mg/L

Ürün Numarası:

2471100

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



The most widely used disinfectant for drinking water, chlorine is also important for sanitizing swimming pools, cooling towers, other industrial equipment, and in the treatment of municipal wastewater. Its measurement and control are vital for both safety and economic reasons. Kit utilizes the Hach Digital Titrator and contains 100 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Kit contains all required reagents and apparatus to perform analysis in the field

Teknik özellikler

Dipnot:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
Metot:	Digital Titrator/Thiosulfate
Muhafaza Türü:	F
Ölçüm aralığı:	20 - 2,000 mg/L Cl ₂
Parametre:	Chlorine, mid and high range - As Total Cl ₂
Platform:	Dijital Titratör
Smallest increment:	0.2 - 5.0
Test kiti modeli:	CN-DT
Test sayısı:	100