



AccuVac kit, dissolved oxygen

Ürün Numarası:

2515050

TRY Fiyat:

Bize Ulaşın

3 gün içinde teslimat

Tehlikeli



Low DO levels usually indicate serious pollution. Adequate amounts are crucial for fish life, but conversely, DO must be excluded from boiler feed water to prevent corrosion. The simple, low-cost color disc method is the most accurate visual test kit method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Teknik özellikler

Metot:	Color disc/HRDO
Model:	Color disc/HRDO
Muhafaza Türü:	G
Ölçüm aralığı:	0 - 15 mg/L O ₂
Parametre:	Dissolved Oxygen - As DO
Platform:	Accuvac
Smallest increment:	0.2
Test kiti modeli:	AccuVac® Kit
Test sayısı:	25