



Sulfite Test Kit, Model SU-5

Ürün Numarası:

148002

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



Utilizing drop count titration/iodometric to determine sulfite, this kit contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Teknik özellikler

Metot:	Drop count, titration Iodometric
Muhafaza Türü:	D
Ölçüm aralığı:	1.0 - 20 mg/L & 10 - 200 mg/L SO_3^{2-}
Parametre:	Sulfite
Platform:	Drop Count Titration
Smallest increment:	1 10
Test kiti modeli:	SU-5
Test sayısı:	100