



Test kit, sulfite, model SU-DT

Ürün Numarası:

2063300

TRY Fiyat:

Bize Ulaşın

2 hafta içinde teslimat

Tehlikeli



Using the Digital Titrator, this kit contains 100 tests with a range of 10-800 mg/L.

Digital Titrator and standardized titration cartridges provide laboratory accuracy for field titrations

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Teknik özellikler

En küçük artış adımları:	0,4-4,0**
Metot:	Digital Titrator/Iodometric
Muhafaza Türü:	F
Ölçüm aralığı:	4 - 800 mg/L SO_3^{2-}
Parametre:	Sulfite
Platform:	Dijital Titratör
Test kiti modeli:	SU-DT
Test sayısı:	100