



Be Right™



Test kit, sulfate, model SF-1

Ürün Numarası:

225100

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



Widely found in natural waters, sulfate may be at high levels in mine drainage. High concentrations of magnesium or sodium sulfate in water act as cathartics. This test kit provides economical determination of sulfate levels using the Extinction/Turbidimetric methods.

Economical test kit for determination of sulfate levels using the Extinction/Turbidimetric Method

Teknik özellikler

Metot:	Extinction/Turbidimetric
Muhafaza Türü:	D
Ölçüm aralığı:	50 - 200 mg/L
Parametre:	Sulfate
Smallest increment:	50
Test kiti modeli:	SF-1
Test sayısı:	100