



Test kit, alkalinity, model AL-DT

Ürün Numarası:

2063700

Tehlikeli

TRY Fiyat:

Bize Ulaşın

2 hafta içinde teslimat



Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Determine phenolphthalein and total alkalinity by visual endpoint titration. Using the Digital Titrator, this contains 100 tests with a range of 10-4,000 mg/L.

Digital Titrator and standardized titration cartridges provide laboratory accuracy for field titrations

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Teknik özellikler

Dipnot:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
En küçük artış adımları:	0.1-10
Metot:	Digital Titrator/Sulfuric acid
Muhafaza Türü:	F
Ölçüm aralığı:	10 - 4,000 mg/L CaCO ₃
Parametre:	Alkalinity - (P) & Total (MO) - As CaCO ₃
Platform:	Dijital Titratör
Test kiti modeli:	AL-DT
Test sayısı:	100