



Bottle, BOD, 300 mL, 6/pk

Ürün Numarası:

62106

TRY Fiyat:
Mevcut

Bize Ulaşın

Set of 6 glass BOD bottles, 300 mL, with stoppers. For use with LBOD sensor to determine the Biological Oxygen Demand. Borosilicate clear glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stoppers included.

Teknik özellikler

Body Material:	Glass
Hacim / Paket boyutu:	6 /pk
Kapasite:	300 mL
Malzeme:	Borosilicate Glass
Metot numarası:	8043
Package size UOM:	6pk