



Beaker, polypropylene, low form, 400 mL, 6/pk

Ürün Numarası: 108078
TRY Fiyat: Bize Ulaşın
Mevcut

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Teknik özellikler

| | |
|------------------|-------------------------|
| Birim: | 6 /pk |
| Çap: | 99.1 mm with spout |
| Gövde malzemesi: | Polypropylene |
| Kapasite: | 400 mL |
| Malzeme: | Polypropylene, Low Form |