



Calcium Standard Solution, 100 mg/L, 500 mL

Ürün Numarası: 4457649
TRY Fiyat: Bize Ulaşın
Mevcut

Ion Selective Electrode (ISE) standard solutions for laboratory analyses

ISE Standard Solution. Calcium, 100 mg/L ± 1 mg/L Ca^{2+} (NIST). 500 mL bottle

Traceable to NIST standard reference materials

Always performing measurement results you can trust is essential. The quality of the standard solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate ISE measurements is good care and calibration of the Ion Selective Electrode.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Teknik özellikler

Hacim / Paket boyutu:	500 mL
Konsantrasyon:	100 mg/L ± 1 mg/L Ca^{2+}
Önemli açıklamalar:	Observed shelf life after opening, following recommended storage and use: 3 months
Package settings:	Bottle
Parametre:	Calcium
Platform:	Çözelti
Saklama koşulları:	10 - 25 °C (protect from light)