



Flask, volumetric, polypropylene, 100 mL

Ürün Numarası: 1406042
TRY Fiyat: Bize Ulaşın
Mevcut

Polypropylene flask and screw cap. Durable with excellent chemical resistance, with the same accuracy as glass. Permanent color graduation line with marking spot. Calibrated “to contain/deliver” at 20°C. The 200/203 capacity flask has dual marks useful in DO determination via Winkler Titration.

Teknik özellikler

| | |
|------------------|----------------|
| Birim: | 1 each |
| Gövde malzemesi: | Polypropylene |
| Kapasite: | 100 mL |
| Sertifikalar: | Not Applicable |