



Funnel, analytical, 65 mL approx. volume

Ürün Numarası:	108367
TRY Fiyat:	Bize Ulaşın
Mevcut	

Standard long stem funnel meets the needs of a wide variety of applications. Mounded ribbing prevents "air lock" commonly occurring during funneling and filtrations. A steep, 58-degree cone slope makes for fast and easy flow of either liquids or small solids. Polypropylene construction is resistant to chemicals and makes this funnel fully autoclavable for sterilization.

Teknik özellikler

Hacim / Paket boyutu:	1 each
-----------------------	--------