



Cylinder, Graduated, 25 mL, Certified

Ürün Numarası: 2636240
TRY Fiyat: Bize Ulaşın
2 hafta içinde teslimat

Certified Class A, this glass cylinder is calibrated "to deliver" according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in white marking. Bumper guard included for all sizes $\geq$ 25 mL.

Teknik özellikler

Gövde malzemesi:	Glass
Hacim / Paket boyutu:	1 each
Kapasite:	25 mL
Ölçüm metodu:	To Deliver
Sertifikalar:	Class A
Tür:	Tek Metrik
Yükseklik:	192 mm