



Beaker, griffin, low form, glass, 50 mL, 12/pk

Ürün Numarası: 50071H
TRY Fiyat: Bize Ulaşın
Mevcut

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologic

Teknik özellikler

Birim:	12 /pk
Gövde malzemesi:	Glass
Kapasite:	50 mL
Malzeme:	Glass, Griffin, Low Form