



Bottle, BOD, black, 300 mL

Ürün Numarası: 2086047
TRY Fiyat: Bize Ulaşın
Teslimat tarihi için arayın

The standard for BOD analysis. Borosilicate glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stopper included (60 and 300 mL sizes). Serialized bottles have silk-screened, sequential

Teknik özellikler

Body Material:	Glass
Hacim / Paket boyutu:	1 each
Kapasite:	300 mL
Malzeme:	Borosilicate Glass
Package size UOM:	1pk