



## Sampler, water core, "sludge judge" large capacity, 1-1/4" diameter

Ürün Numarası:

2506102

TRY Fiyat:  
Mevcut

Bize Ulaşın

Accurate readings can be obtained on suspended solids up to a depth of 4.5 m (15'). Comes in three 1.5 m (5') sections with screw-type fittings. Designed for sampling solids in wastewater and food processing applications. Standard Sludge Judge features 3/4" diameter plastic tubing, while the large capacity Sludge Judge features 1-1/4" diameter plastic tubing for double the capacity.

### Teknik özellikler

Çap:	31.75 mm
Hacim / Paket boyutu:	1 each