



Desiccator, vacuum, 171 mm (6.7") diameter

Ürün Numarası: 2238300
TRY Fiyat: Bize Ulaşın
Teslimat tarihi için arayın

Round polycarbonate top provides lightweight vessel with unmatched transparency. Polypropylene lower section creates support chamber at full vacuum. Stopcock for vacuum applications (1/4-3/8" OD vacuum tubing, sold separately). Holds 28 in Hg (711 mm) for 24 hrs. Silicone O-ring makes grease unnecessary. For room temperature, non-organic solvent, and non-autoclave uses only.

Teknik özellikler

Hacim / Paket boyutu: 1 each