



Pipet rack

Ürün Numarası:

2094700

TRY Fiyat:
Mevcut

Bize Ulaşın

Polypropylene, fifty-place rack for safe, orderly storage of pipets. Helps create an efficient storage system or a well organized liquid handling station. Holes in top and middle decks are 16 mm diameter (0.63"); bottom holes are tapered to 6 mm (0.24") to prevent contact with supporting surface. Autoclavable to 122°C. Overall dimension: 21 x 11 x 22 cm (8.3 x 4.3 x 8.7").

Teknik özellikler

Derinlik:	8.7 Inches
Genişlik:	4.3 Inches
Hacim / Paket boyutu:	1 each
Uzunluk:	210 mm