



Test Tube, rack, polyethylene, 20 mm (O.D.), 20 holes

Ürün Numarası:	2497912
TRY Fiyat:	Bize Ulaşın
Mevcut	

Twenty-place rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Autoclavable. Stackable and autoclavable. Overall dimensions: 24.6 x 10.7 x 6.4 cm (9.7 x 4.2 x 2.5"). Color may Vary.

Teknik özellikler

Çap:	20 mm
Hacim / Paket boyutu:	1 each
Test tüpü düzenlemesi:	4 x 5
Test tüpü yuvası:	20