



## Iron and manganese color disc test kit, model IR-20

Ürün Numarası:

146300

TRY Fiyat:  
Mevcut

Bize Ulaşın

Tehlikeli



When customers complain of stained and dingy laundry, this kit lets you do quick, on-the-spot tests for iron and manganese to see where the problem lies.

Method / Chemistry: Color disc colorimetry. Test: Iron/ Manganese. Range (mg/L): 0-5/ 0-3.

**Simple, fast, color disc method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

### Teknik özellikler

En küçük artış adımları:	0.2
Model:	IR-20
Ölçüm aralığı:	0 - 3,0
Ölçüm aralığı (2):	0 - 4,0 mg/L Fe
Parametre:	Iron and Manganese
Platform:	Renk Diski
Uygulama:	Water Quality