



## Test kit, carbon dioxide, model CA-DT

Ürün Numarası:

2064100

Tehlikeli

TRY Fiyat:

Bize Ulaşın

Teslimat tarihi için arayın



Utilizing the Hach Digital Titrator, the Model CA-DT kit is used to determine carbon dioxide concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10.1-1,000 mg/L. Contains approximately 50 tests.

**Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration**

**Rugged titration device is ideal for field use, with no glass burets required**

### Teknik özellikler

Dipnot:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
Metot:	Digital Titrator/  Sodium hydroxide  Phenolphthalein
Muhafaza Türü:	F
Ölçüm aralığı:	10 - 1,000 mg/L CO <sub>2</sub>
Parametre:	Carbon Dioxide - As CO <sub>2</sub>
Platform:	Dijital Titratör
Smallest increment:	0.1, 2.0
Test kiti modeli:	CA-DT
Test sayısı:	50