



Test kit, acidity, model AC-DT

Ürün Numarası:

2064000

TRY Fiyat:

Bize Ulaşın

2 hafta içinde teslimat

Tehlikeli



Acidity is usually caused by the presence of mineral acids, salts of strong acids, or free carbon dioxide. In surface or ground waters, acidity may form natural substances or industrial pollution. Excessive levels cause corrosion and can be detrimental to fish. Utilizing the Hach Digital Titrator, the Model AC-DT kit is used to determine acidity concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10-4,000 mg/L. Contains approximately 100 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Teknik özellikler

Dipnot:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
En küçük artış adımları:	0.1-1.10
Metot:	Digital Titrator/Sodium hydroxide
Muhafaza Türü:	F
Ölçüm aralığı:	10 - 4,000 mg/L CaCO ₃
Parametre:	Acidity - (P or Total) - As CaCO ₃
Platform:	Dijital Titratör
Test kiti modeli:	AC-DT
Test sayısı:	100