



Test kit, hypochlorite, model CN-HRDT

Ürün Numarası:

2687100

Tehlikeli

TRY Fiyat:

Bize Ulaşın

3 gün içinde teslimat



Hypochlorite Kit, Model CN-HRDT with Digital Titrator

Utilizing the Hach Digital Titrator, this kit is used to determine hypochlorite, high range (Bleach) solution concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 5 - 15% as Cl₂: 50 - 150 g/L as Cl₂. Contains 100 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Designed for on-site analysis of bulk hypochlorite deliveries or bulk storage solutions at drinking water or chloramination facilities

Verify concentration of hypochlorite solutions at time of delivery

Kit contains all required reagents, solutions, and glassware

Teknik özellikler

Metot:	Digital Titrator/ Thiosulfate
Metot numarası:	10100
Muhafaza Türü:	G
Ölçüm aralığı:	5 - 15 % = 50-150 g/L Cl ₂
Parametre:	Hypochlorite, high range (Bleach) - As Cl ₂
Platform:	Dijital Titrator
Smallest increment:	0.05%
Test kiti modeli:	CN-HRDT
Test sayısı:	100