



Be Right™



Bottle, sampling, polycarbonate, 500 mL, 4/pk

Ürün Numarası:	2324539
TRY Fiyat:	Bize Ulaşın
Mevcut	

Transparent polycarbonate bottle with cap features wide mouth (3.8cm/1.48" diameter) opening for easy filling and square design for maximum space efficiency. Autoclavable with molded-in graduations (do not engage threads for autoclave cycle.)

Teknik özellikler

Body Material:	Polycarbonate
Hacim / Paket boyutu:	4 /pk
Kapasite:	500 mL