



Filter funnel, micro

Ürün Numarası:

2619000

TRY Fiyat:
Mevcut

Bize Ulaşın

Feature a completely assembled filtration unit: Funnel, filter, and petri dish. Filter: 47 mm, white with gridded surface, 0.45 µm GN-6 Metrical media Funnel: Polypropylene Petri Dish: polystyrene lid and polypropylene base Fits 28315-00, 28461-00 and 28461-50 manifolds. Sterile and individually bagged.

Teknik özellikler

Çap:	47 mm
Gözenek boyutu:	0.45 µm
Hacim / Paket boyutu:	50
Malzeme:	Mixed-Cellulose Ester