



Be Right™



Filter, membrane, pore size 0.2 µm, diameter 47 mm, 200/pk (Gelman)

Ürün Numarası: 2618800
TRY Fiyat: Bize Ulaşın
2 hafta içinde teslimat

For isolation and enumeration of stressed organisms (Pseudomonas species). High retention and flow rate. Certified for use in membrane filtration analysis. White polyethersulfone media. STerile and individually packaged. Pore size: 0.2 µm. Diameter: 47mm.

Teknik özellikler

Çap:	47 mm
Filter Type:	Filter
Gözenek boyutu:	0.2 µm
Hacim / Paket boyutu:	200 /pk
Malzeme:	Polyethersulfone