



Sample vial, 120 mL w/out thiosulfate, 200/pk

Ürün Numarası: 2634002
TRY Fiyat: Bize Ulaşın
Mevcut

120 mL sterile, hinged-cap vials meet all EPA requirements for microbiology sampling. Collection, incubation, and accurate reading can all be done in the vial without a transfer step. 100 mL raised fill line. Polypropylene.

Teknik özellikler

Çap: 43 mm
Hacim / Paket boyutu: 200/pk
Kapasite: Fill line at 100 mL ($\pm$ 2.5 mL)
Önemli özellik: Sample Vial, without Thiosulfate
Uzunluk: 93 mm