



## Acid solution for EZ4003 and EZ4004 Total Alkalinity Analyser, 0.025 N HCl (10-250 mg/L), 2 L

Ürün Numarası:

APPC4004-01-EU

TRY Fiyat:

Bize Ulaşın

### Ready-to-use Reagents for EZ Series Online Analysers

#### Ready-to-use

Hach EZ reagents are ready-to use with EZ Series Online Analysers. Since there is no additional preparation needed this eliminates an additional error source when using the EZ instruments and reagents system.

#### Confidence in your results

Reduce the sources of errors in your workflow to increase the confidence in your results.

#### Less complexity

Using ready-to-use EZ reagent kits from Hach can reduce the number of handling steps around your analytical equipment significantly which allows less hands-on-time and complexity.

---

### Teknik özellikler

Hacim / Paket boyutu:

2 L

Konsantrasyon:

0.025 N HCl

Parametre:

Total Alkalinity