



pH Tampon Çözelti, pH 7,00 (500 mL)

Ürün Numarası: 2283549
TRY Fiyat: Bize Ulaşın
Mevcut

pH Buffer standard solutions for laboratory and field analyses

Kaliteli Tampon Çözelti, kolay tanımlama için sarı renk kodlu. NIST izlenebilir. pH 7,00 ± 0,02; 25°C'de.

Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard pH Buffer Solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Teknik özellikler

Hacim / Paket boyutu:	500 mL
Konsantrasyon:	7,00 ± 0,02 pH
NIST izlenebilir:	Evet
Önemli açıklamalar:	Observed shelf-life after opening, following recommended storage and use: 6 months
Package settings:	Bottle
pH Değer(ler)i:	7.00 ±0.02 pH at 25°C
Renk kodlu?:	Evet
Saklama koşulları:	10 - 25 °C (protect from light)
Tanım:	Tampon Çözelti, (Sarı renk kodlu) (NIST)