



Free and total chlorine test strips, 0 - 10 mg/L, 50 tests

Ürün Numarası: 2745050
TRY Fiyat: Bize Ulaşın
Mevcut

In the field, technicians can test many samples in only a few minutes, and make immediate evaluations on-site. No measuring, set-up, clean-up, or chemical handling are necessary. Lab technicians also use the test strips for pre-test screening-to detect the presence of materials that might interfere with testing. Steps: (mg/L) 0, 0.5, 1.0, 2.0, 4.0, 10.0.

Teknik özellikler

Birim:	50 tests
Dipnot:	*mg/L unless otherwise noted; ppb = µg/L; ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
En küçük artış adımları:	Steps: 0, 0.5, 1, 2, 4, 10
Ölçüm aralığı:	0 - 10 mg/L Cl ₂
Parametre:	Free and Total Chlorine
Platform:	Test Strip
Test kiti modeli:	TS