



Test strips pH , 4-9 pH units

Ürün Numarası:	2745650
TRY Fiyat:	Bize Ulaşın
Mevcut	

Most natural waters range from pH 4 to pH 9, but commonly are above pH 7 because of carbonates and bicarbonates (alkalinity). The pH of water used in industry, boiler feed water, and swimming pools usually is kept within a narrow range. Contains 50 tests.

Easy to use, disposable, and inexpensive

Great for quick evaluations of pH levels in the lab or in the field

Teknik özellikler

Birim:	50 tests
En küçük artış adımları:	Steps: 4, 5, 6, 7, 8, 9 pH
Ölçüm aralığı:	4 - 9 pH units
Parametre:	pH, mid range
Platform:	Test Strip
Test kiti modeli:	Test Strips