



Desiccator cabinet, non-vacuum

Ürün Numarası:

2618100

TRY Fiyat:
Mevcut

Bize Ulaşın

Stackable, polystyrene, transparent from any side. Door has foam rubber gasket, two roller latches for airtight fit. Inside dimensions: 21 x 18 x 13 cm (8.25 x 7 x 5"). Door opening: 19 x 11 cm (7.5 x 4.4"). Desiccator Grease and silica gel available and sold separately.

Teknik özellikler

Hacim / Paket boyutu:

1 each