



Soxhlet thimbles, cellulose, 33 x 80 mm, 25 pk

Ürün Numarası: 2966700

ÜRETİMDEN KALDIRILAN ÜRÜN
Bu ürün artık mevcut değil.

Xenosep extraction thimbles are ideal for determining the concentration of nonvolatile and semivolatile organic compounds from environmental solid waste matrices, sludge, soil, air, and water.

Xenosep extraction thimbles are ideal for determining the concentration of nonvolatile and semivolatile organic compounds from environmental solid waste matrices, sludge, soil, air, and water.

The thimbles are manufactured from high purity alpha cellulose cotton linters which provide excellent mechanical strength, stability, and retention with very low/no background extractables in a disposable, sample holder format.

Teknik özellikler

Gözenek boyutu:	10 µm
Hacim / Paket boyutu:	25 /pk
Malzeme:	Cellulose