



Tube, culture, w/ screw cap, 25 x 150 mm, 10/pk

Ürün Numarası: 2096250
TRY Fiyat: Bize Ulaşın
Mevcut

The cap has a firmly bonded rubber liner, providing a sure and highly inert sealing surface.

Cap has a firmly bonded rubber liner, providing a sure and highly inert sealing surface.

Teknik özellikler

Boyutlar (D x U):	25 mm x 150 mm
Çap:	25 mm
Hacim / Paket boyutu:	10 /pk
Kapasite:	Capacity to Neck 50 mL
Önemli özellik:	Tube, 25 x 150 mm
Uzunluk:	150 mm