



Thermometer, 30 to 50 C°, 250 mm

Ürün Numarası:

2999200

TRY Fiyat:

Bize Ulaşın

Mevcut

Non-mercury thermometer, 30 - 50 degrees C Length: 250 mm Immersion: 76 mm Safety coated.

Blue fill solution on white background provides excellent visibility.

Non-mercury thermometer, 30 - 50 degrees C Length: 250 mm Immersion: 76 mm Safety coated.

Blue fill solution on white background provides excellent visibility.

Teknik özellikler

Hacim / Paket boyutu:	1
Sıcaklık aralığı:	30 °C - 50 °C
Termometre Akışkan Türü:	Non-Mercury, blue
Termometre taksimatları:	0.1 °C
Uzunluk:	250 mm