



Bottle, sampling, presterilized, round, 50/pk

Ürün Numarası: 2495050
TRY Fiyat: Bize Ulaşın
Mevcut

Polystyrene, disposable, transparent, non-fluorescent, break-resistant, presterilized bottles with caps. Bottles have 100 mL fill-to line; capacity of 175 mL. Bottles used with medium containing MUG or other fluorescing reagent. Work well with P-A Broth Ampules and PathoScreen Medium.

Teknik özellikler

Dechlorinating Agent?: No
Hacim / Paket boyutu: 50/pk