



SAC950 Sample changer for automation of large sample series (Radiometer Analytical)

Ürün Numarası: R41T044

ÜRETİMDEN KALDIRILAN ÜRÜN

Bu ürün artık mevcut değil.

Effortless sample handling due to total automation

Extended sample changer for automation of large sample series, multi-radius sample trays for up to 70 samples

Compact and adaptable

Flexible and programmable

Quick and efficient

Safe and clean

Teknik özellikler

Ağırlık:	27 kg
Boyut (Y x G x D):	70 mm x 580 mm x 686 mm
Tür:	Numune Değiştirici
Uzunluk:	686 mm